

WIC Breast Pump Guide:

| | |
|---|---|
|  | <p>MANUAL OR HAND PUMP: (http://www.medelabreastfeedingus.com/products/breast-pumps/157/harmony-breastpump)</p> <ul style="list-style-type: none">■ For occasional or short term pumping (i.e. stay at home mothers, mild engorgement, short separations, etc.).■ Unit contains a standard 24 mm shield (flange) and a 27 mm shield option.■ Not intended for reuse, not to be returned. |
|  | <p>PERSONAL DOUBLE ELECTRIC PUMP (AMEDA) : (http://www.ameda.com/products/personal.aspx)</p> <ul style="list-style-type: none">■ For mothers of premature infants or infants with serious or chronic medical problems and mothers returning to work or school.■ Unit contains a standard 24 mm shield (flange) and a 27 mm shield option.■ Not intended for reuse, not to be returned. |
|  | <p>HOSPITAL GRADE ELECTRONIC PUMP (LACTINA): (http://www.medelabreastfeedingus.com/products/breast-pumps/74/lactina-select-breastpump)</p> <ul style="list-style-type: none">■ Loan for short-term medical complications of mother or infant (i.e. sever engorgement, nipple trauma, temporary latch-on problem, surgery, separation, etc.).■ Reusable loaner pump; property of LHU. Must be cleaned with 1:10 bleach solution after each loan. |